



How to Take Care of Your Compression Stockings

Your new compression stockings should last approximately one year with proper care. Wash daily by hand, or in a laundry bag in a washing machine (delicate wash cycle) on cold or lukewarm. We recommend using a mild detergent such as a baby detergent or SIGVARIS washing solution. Press the stockings between two dry towels and lay them out to dry in a shady place with plenty of ventilation.

Do not use fabric softener.

Do not put the stockings in the dryer as this will cause them to overheat break down the fibers.

Do not dry clean.

Do not use chemicals.

Do not place on heaters/radiators.

Do not iron

Do not cut off any loose threads inside the stockings.

Sharp fingernails and rings with sharp edges can quickly result in defects. Pay attention to these risks when donning and removing your stockings or use SIGVARIS donning gloves.